

Multi-Purpose Teaching Flume - C4-MKII



The Armfield Multi-purpose Teaching Flume has been specifically designed to demonstrate the principles of fluid mechanics when applied to engineering structures in open channel flow.

The MKII version is supplied as an accessory to the F1-10 Hydraulics Bench and can be used with any recently supplied F1-10.

COST EFFECTIVE LABORATORY FLUME FOR TEACHING AND RESEARCH PURPOSES



C4-MKII - Teaching flume



C4-69 - Artificially roughened bed



Sharp Crested Weir (in vented condition)

Experimental content

- ▶ Study of open channel flume
- ▶ Use of hook and point gauges to measure water level
- ▶ Learning how to apply force-momentum and steady-flow energy equations to simple flow situations
- ▶ Understanding the relationship between water level above the crest of a weir and flow rate over the weir
- ▶ Using hydraulic structures to control level, eg syphon spillways
- ▶ Understanding sub and super-critical flow and the underlying characteristics of waves
- ▶ Hydraulic jump
- ▶ Using hydraulic structures for control of flow, eg sluice gate
- ▶ Applying and understanding Manning's formula
- ▶ Use of a Pitot static tube to measure flow rate and velocity profiles (using optional C4-61)

UK office - email: sales@armfield.co.uk tel: +44 (0) 1425 478781 (for ROW)
 USA office - email: info@armfield.inc tel: +1 (609) 208-2800 (USA only)

Issue: 2
 URL: <http://www.armfield.co.uk/c4mkii>

Applications
ChE ME CE IP

We reserve the right to amend these specifications without prior notice. E&OE © 2019 Armfield Ltd. All Rights Reserved

Description

The C4MKII is a small open channel flume, available in 2.5m or 5.0m lengths, with clear acrylic sides to the working section for complete visibility of the flow.

The channel is fitted with a PVC inlet tank and is designed for free discharge into the Hydraulics Bench. The flume is mounted on a rigid framework and can be tilted by use of a calibrated screwjack, which enables accurate slope adjustment of the channel.

The inlet tank incorporates a stilling arrangement to diffuse the water flow prior to entry into the channel, ensuring smooth uniform flow. The level in the working section of the flume is controlled using an overshoot weir (stop logs) at the discharge end.

Bed-pressure tappings and fixing points for models are provided.

A longitudinal scale positioned at the top of the channel enables depth gauges and pitot static tubes to be accurately positioned along the channel length.

The flume is designed for use with a standard Armfield F1-10 Hydraulics Bench, which provides the pumped water flow, the flow control valve and a volumetric tank for flow measurement.


Also available is an optional flow meter, which can be fitted to the C4-MKII to enable direct flow measurements to be taken.

Optional educational software is available (C4-MKII-ABASIC) offering a complete teaching package of coursework. The student manually enters data in the software, which can then be used for calculations, data processing and graph plotting.

Detailed experimental content

- ▶ To determine the relationship between upstream head and flowrate for water flowing over a Sharp Crested weir
- ▶ To calculate the discharge coefficient and to observe the flow patterns obtained over a Sharp Crested weir
- ▶ To determine the relationship between upstream head and flowrate for water flowing over a Broad Crested weir (long base weir)
- ▶ To calculate the discharge coefficient and to observe the flow patterns obtained over a Broad Crested weir (long base weir)
- ▶ To determine the relationship between upstream head and flowrate for water flowing over a Crump weir
- ▶ To determine the modular limit and to observe the flow patterns obtained over a Crump weir
- ▶ To determine the relationship between upstream head and flowrate for water flowing under a sluice gate (undershot weir)
- ▶ To calculate the discharge coefficient and to observe the flow patterns obtained under a sluice gate (undershot weir)
- ▶ To determine the relationship between upstream head and thrust on a sluice gate (undershot weir) for water flowing under the sluice gate
- ▶ To determine the relationship between the specific energy and upstream head for water flowing under an undershot weir
- ▶ To investigate the characteristics of a standing wave (the hydraulic jump) produced when water flows beneath an undershot weir and to observe the flow patterns obtained
- ▶ To determine the relationship between upstream head and flowrate for water flowing through a Venturi flume
- ▶ To calculate the discharge coefficient and to observe the flow patterns obtained through a Venturi flume
- ▶ To determine the characteristics and observe the flow patterns obtained for water flowing through a culvert. Requires C4-62
- ▶ To observe the flow patterns obtained for water flowing around splitters with different profiles. Requires C4-63
- ▶ To observe the flow patterns associated with the flow of water over a dam spillway. Requires C4-64
- ▶ To determine the relationship between upstream head and flowrate through a siphon spillway in the "blackwater" fully primed condition. Requires C4-65
- ▶ To calculate the discharge coefficient and to observe the operation of the siphon as it primes and de-primed. Requires C4-65
- ▶ To determine the relationship between upstream head and flowrate through an air regulated siphon. Requires C4-65
- ▶ To calculate the discharge coefficient and to observe the operation of the air regulated siphon as it primes and de-primed. Requires C4-65
- ▶ To determine the relationship between upstream head and flowrate beneath a radial gate (Tainter Gate) under different operating conditions. Requires C4-66
- ▶ Calculate the discharge coefficient in each position of the radial gate. Requires C4-66
- ▶ To observe the flow patterns associated with the flow of water over different bed profiles. Requires C4-68
- ▶ To determine the effect of a roughened bed on the depth of water at different flowrates and to obtain appropriate coefficients to satisfy the Manning Formula. Requires C4-69
- ▶ To determine the effect of different wave regimes on different flow patterns. Requires C4-67

Accessories supplied with the standard flume



Other available accessories for the flume


C4-65 - Air-regulated syphon



C4-64 - Free overflow spillway section complete with ski jump, sloping apron and blended reverse curvature attachments



C4-65-Syphon spillway



C4-61 - Pitot tube and manometer



C4-63 - Flow splitters; central wall with various nose pieces




C4-66 - Model radial gate



C4-68 - False floor sections



C4-62 - Culvert fitting has square, rounded edges



C4-67 - Wave generator and wave absorbing beach




C4-69 - Artificially roughened bed



Requirements

Scale

F1-10 PC USB



- ▶ Electrical supply: The F1-10 requires an electrical supply, please refer to the F1 data sheet for details.
- ▶ F1-10-A: 220-240V / 1ph / 50Hz / 10A
- ▶ F1-10-B: 110-120V / 1ph / 60Hz / 10A
- ▶ F1-10-G: 220V / 1ph / 60Hz / 10A
- ▶ F1-10-2-A: 220-240V / 1ph / 50Hz / 10A
- ▶ F1-10-2-B: 110-120V / 1ph / 60Hz / 10A
- ▶ F1-10-2-G: 220V / 1ph / 60Hz / 10A



Technical specifications

Models and gauges supplied

Venturi flume
Sharp and broad-crested weirs
Crump weir
Adjustable undershot weir
Two Vernier level gauges (Hook and point gauges)

Channel dimensions

Width	76mm
Height	250mm
Channel Slope	Adjustable between -1% and +3%

Overall dimensions

Model	C4-MKII-2.5M	C4-MKII-5.0M
Length	2.91m	5.41m
Width	0.62m	0.62m
Height	1.46m	1.46m

Packed and crated shipping specifications

Volume	2.8m ³	3.0m ³
Gross weight	330kg	380kg

Related products

- C4-61: Pitot tube and manometer
- C4-62: Culvert fitting, one edge square, one rounded
- C4-63: Flow splitters; central wall with various nose pieces
- C4-64: Free overflow spillway section complete with ski jump, sloping apron and blended reverse curvature attachments
- C4-65: Syphon spillway and air-regulated syphon
- C4-66: Model radial gate
- C4-67: Wave generator and wave absorbing beach
- C4-68: False floor sections for gradually varied profiles
- C4-69: Artificially roughened bed 2.5m-long section (two required for a 5m flume)

Essential accessories / equipment

Armfield F1-10 Hydraulics Bench

Software

Optional educational software is available (C4-MKII-ABASIC) offering a complete teaching package of coursework. The student manually enters data in the software, which can then be used for calculations, data processing and graph plotting.

Ordering specification

- ▶ A 76mm wide, 250mm high open channel for use with an F1-10 Hydraulics Bench
- ▶ Available in 2.5m and 5.0m working section lengths
- ▶ Clear acrylic sides to give visibility of the working section
- ▶ A jacking system permits the slope of the channel bed to be adjusted between -1% and +3%
- ▶ Inlet tank with flow stilling arrangement
- ▶ Includes a Venturi, sharp and broad-crested weirs, 2 vernier level gauges, adjustable undershot weir and crump weir
- ▶ Wide range of other models available as accessories.
- ▶ Optional flow meter
- ▶ Comprehensive instruction manual

Ordering codes

- C4-MKII-2.5M-10**
Multi-purpose Flume, 2.5m, with basic accessories
- C4-MKII-5.0M-10**
Multi-purpose Flume, 5m, with basic accessories
- C4-MKII-2.5M-11**
Multi-purpose Flume, 2.5m, with basic accessories and flow meter
- C4-MKII-5.0M-11**
Multi-purpose Flume, 5m, with basic accessories and flow meter
- C4-MKII-ABASIC**
Educational Software

Armfield standard warranty applies with this product

Knowledge base

- > 28 years expertise in research & development technology
- > 50 years providing engaging engineering teaching equipment

Benefit from our experience, just call or email to discuss your laboratory needs, latest project or application.

An ISO 9001:2015 Company



armfield.co.uk

Aftercare

Installation
Commissioning
Training
Service and maintenance
Support: armfieldassist.com